




USB Heating Eye Mask iFlo VEM-200

Provide Effective Warm Compress to Dry Eyes

Product Concept

iFlo is a smart eye mask designed to provide heating comfort to the eyes. It works by connecting eye mask with power bank or computer or any other power source. It can relieve eye fatigue and discomfort caused by dry eye or other eye problems. It can also accelerate the blood circulation around the eye to improve black eye circle conditions. There is a removable sheet on the eye mask for users to clean and keep hygiene. The sheet is machine washable and dryable.



Temperature
Adjustable



Nasal Arch Fit



Machine Washable
Sheet




Auto Switch off



FlexStrip™ Fits Face Better

FlexStrip™ is innovatively designed and applied in the nose bridge area of the iFlo eye mask. FlexStrip™ is adjustable according to the nose shape to make the iFlo eye mask better fit on the face. Better fit provides better heating effect at the eye lid area.



- 44°C/111°F*
- 41°C/106°F
- 39°C/102°F


3 Heating Options with Smart Control

iFlo provides 3 optional heating temperatures to meet different comfort demands. The heating time is 12 minutes for each option. When the 12 minutes' heating SPA is finished, it will automatically switch off.

*The temperature varies at $\pm 3^{\circ}\text{C}$ / $\pm 37^{\circ}\text{F}$ depending on the environmental temperature.

High Quality Material in Each Layer

Each layer of the eye mask is made up with high quality materials to ensure comfortable and safe wearing experience.



Specification

VEM-200

Heating method	Powered by USB connection
Heating temperature	39°C/41°C/44°C, ±3°C (102°F/106°F/111°F, ±3°F)
Heating duration	12 minutes
Life span	750 Heating times (Max)
Input voltage	DC 5V/1A
Power	4W (max.)
Color	Space grey
size	22cm x 7cm
Weight	25g

VisuScience Meditech Co., Ltd.

Add: Room 718, No. 208, Guangyi Road, Wuxi Jiangsu 214016, China

Tel: +86-510-85757880

E-mail: info@visuscience.com

Web: www.visuscience.com



Scan to Learn More