



USB Heating Eye Mask iFlo VEM-200

Provide Effective Warm Compress to Dry Eyes

Listen to the clinical

iFlo VEM-200

USB Heating Eye Mask

Provide Effective Warm Compress to Dry Eyes

- Fit face better with adjustable strip in the middle
- Switch off automatically after 12 minutes heating
- Removable sheet can be washed by machine



○ 39°C/102°F ○ 41°C/106°F ○ 44°C/111°F



VisuScience Meditech Co., Ltd.

Add: Room 501, Ruiyun 5 No. 99 Furongzhongsan Rd. Wuxi 214192, Jiangsu Province, China

Tel: +86-510-85757880 | E-mail: info@visu-science.com | Web: www.visu-science.com



Scan for Details